

## One-stop PURCHASE


Perfect price-performance ratio products

## Professional SER VICE



GD 360 Automatic Preparation System for Test Specimens of Architectural Coatings is suitable for those test specimens applicated by rods. It has the following features:

- ◆ Fully automatic batch preparation of test specimens, almost without manual participation, to meet the needs of large-scale preparation of specimens with high consistency.
- ◆ Specimens digital management is adopted, which can encode and record coating specimens and prepared specimens, so as to realize the traceability management of specimens.
- ♦ According to the requirements of the different standards, operator can finish specimens preparation on Chart, PVC Film and Asbestos-free Cement Panel automatically.
- ◆ Use robot-arm to applicate coating films. This system includes a vacuum drawdown plate and a common plate for cement panels, which can achieve uniform coating film easily on different substrates.
- ♦ Blank test panels is automatically fed according to the demands, and automatically unloaded to the standard specimens tray after finishing coating, which is convenient for maintenance.
- ◆ Automatic cleaning device for applicating tools: after finishing one coating film, the applicating tools will be automatically cleaned, and then the next coating film will be carried out.
- ◆ The system adopts friendly human-computer interface, with intuitive display, convenient operation and real-time monitoring of running status.



## Main Technical Parameters:

- ★ Specification and Capacity of Test Specimens:
  - 250mm × 140mm Charts (50pcs);
  - 150mm × 70mm Asbestos-free Cement Panels (500 pcs)
  - 430mm × 150mm Asbestos-free Cement Panels (100 pcs);
  - $432\text{mm} \times 165\text{mm}$  Black Plastic Panels ( 100 pcs )
- ★ Specimens Preparation Rate: 20 samples per 12 hours (about 400pcs test specimens) or 10 samples per 10 hours (about 200pcs test specimens)
- ★ Single Working Cycle: 10~12 hours (10~20 samples)
- ★ Applicating Tools Capacity: 4 sets of specimens preparation
- ★ Number of Specimens Trays: 20 sets of specimenss
- ★ Cover Area: approximate 30m²
- ★ Total Power: 5 kW
- ★ Ordering Information: BGD 360---Automatic Preparation System for Test Specimens of Architectural Coatings

Tel: +86-20-82169666 E-mail: sale: carish@biuged.com Service: service@biuged.com WhatsApp: +86-13316053991 Website: www.biuged.com