

## Accessories of Viscometer

### THD Series----Low Temperature Thermostatic Bath

#### BGD 1600----Standard Oil (500ml)

#### BGD 1601----Low Viscosity Adapter (The 0th rotor and can measure 1.0 mPa.s sample)

**BGD 1602----Small Sample Adapter** (It's specially designed for customers with small samples and requires only a quantity as small as 5~18ml; Sample cup is easily cleaned, installed, removed, and observed; Accurate data of shear rate and shear stress can be obtained for a minute analysis of the fluid characteristics of the sample; One-off sample cups are available)

#### BGD 1603----Mini Single Color Printer

#### BGD 1604----RTD Thermometer (0°C~100°C)

**BGD 1605----Thermoses** (It has been designed to accurately measure the viscosity of heated oil, paraffin, asphalt emulsions, medicine, high polymer and similar liquid materials at high temperatures. It allows control of a sample's temperature at a range of room temperature +300°C. It can program the intelligent thermostat and ensures its temperature fluctuation within 0.1°C. Sample dosage is as small as 10ml, and one-off sample cups are available. Standard type of SC4 spindle can be applied, which is easy to use)

**BGD 1607----BGD 155 Viscometer data collection and graphing software** (Automated data collection; resulting in viscometer graph, and recording measuring data each time; Allowing comparison of at most 10 history data; Output Excel documents)

**BGD 1608----BGD 156-157-158 Viscometer** data collection and programmed analysis software (Automated data collection and testing; resulting in viscometer graph, and recording measuring data each time; Allowing comparison of at most 10 history data; Output Excel documents; Programmable downloading, Offline working; Analyses data, draw charts and print)




BGD 1601



BGD 1602



BGD 1603



BGD 1607